

THRSP 抗原（重组蛋白）

中文名称：THRSP 抗原（重组蛋白）

英文名称：THRSP Antigen (Recombinant Protein)

别名：S14; Lpgp; THRP; LPGP1; SPOT14

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	thyroid hormone responsive
Synonyms:	S14; Lpgp; THRP; LPGP1; SPOT14
Swissprot:	Q92748
Gene Accession:	BC031989
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is similar to the gene product of S14, a rat gene whose expression is limited to liver and adipose tissue and is controlled by nutritional and hormonal factors. This gene has been shown to be expressed in liver and adipocytes, particularly in lipomatous modules. It is also found to be expressed in lipogenic breast cancers, which suggests a role in controlling tumor lipid metabolism.