

## THYN1 抗原（重组蛋白）

中文名称：THYN1 抗原（重组蛋白）

英文名称： THYN1 Antigen (Recombinant Protein)

别名： MY105; THY28; MDS012; HSPC144; THY28KD

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	thymocyte nuclear protein 1
<b>Synonyms:</b>	MY105; THY28; MDS012; HSPC144; THY28KD
<b>Swissprot:</b>	Q9P016
<b>Gene Accession:</b>	BC006978
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described.