

## TIMM17A 抗原（重组蛋白）

中文名称：TIMM17A 抗原（重组蛋白）

英文名称：TIMM17A Antigen (Recombinant Protein)

别名：TIM17; TIM17A

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	translocase of inner mitochondrial membrane 17A
<b>Synonyms:</b>	TIM17; TIM17A
<b>Swissprot:</b>	Q99595
<b>Gene Accession:</b>	BC000294
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Essential component of the TIM23 complex, a complex that mediates the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.