

兔抗 CA10 多克隆抗体

中文名称：兔抗 CA10 多克隆抗体

英文名称： Anti-CA10 rabbit polyclonal antibody

别名： CARPX; CA-RPX; HUCEP-15

抗原： CA10

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that belongs to the carbonic anhydrase family of zinc metalloenzymes, which catalyze the reversible hydration of carbon dioxide in various biological processes. The protein encoded by this gene is an acatalytic member of the alpha-carbonic anhydrase subgroup, and it is thought to play a role in the central nervous system, especially in brain development. Multiple transcript variants encoding the same protein have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	CA10

Immunogen:	Fusion protein of human CA10
Full name:	carbonic anhydrase X
Synonyms :	CARPX; CA-RPX; HUCEP-15
SwissProt:	Q9NS85
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100

