

TNFAIP8 抗原（重组蛋白）

中文名称：TNFAIP8 抗原（重组蛋白）

英文名称：TNFAIP8 Antigen (Recombinant Protein)

别名：NDED; GG2-1; SCCS2; SCC-S2; MDC-3.13

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	tumor necrosis factor, alpha-induced protein 8
Synonyms:	NDED; GG2-1; SCCS2; SCC-S2; MDC-3.13
Swissprot:	O95379
Gene Accession:	BC007014
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Acts as a negative mediator of apoptosis and may play a role in tumor progression. Suppresses the TNF-mediated apoptosis by inhibiting caspase-8 activity but not the processing of procaspase-8, subsequently resulting in inhibition of BID cleavage and caspase-3 activation.