

兔抗 CAB39L 多克隆抗体

中文名称：兔抗 CAB39L 多克隆抗体

英文名称： Anti-CAB39L rabbit polyclonal antibody

别名： MO2L; MO25-BETA; bA103J18.3; RP11-103J18.3

相关类别： 一抗

抗原： CAB39L

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Calcium-binding protein 39-like is a protein that in humans is encoded by the CAB39L gene. CAB39L (calcium binding protein 39-like), also known as MO25L (MO25-like) or MO2L, is a 337 amino acid protein that is similar to MO25 and is found in the serum of nearly half of all patients diagnosed with acute monocytic leukemia. This suggests a role for CAB39L in carcinogenesis.
Applications:	ELISA, WB, IHC
Name of antibody:	CAB39L
Immunogen:	Fusion protein of human CAB39L
Full name:	calcium binding protein 39-like

Synonyms :	MO2L; MO25-BETA; bA103J18.3; RP11-103J18.3
SwissProt:	Q9H9S4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	39 kDa
WB Positive control:	Human fetal brain tissue and 231 cells
WB Recommended dilution:	500-2000



