

## 兔抗 DDX3Y 多克隆抗体

中文名称：兔抗 DDX3Y 多克隆抗体

英文名称： Anti-DDX3Y rabbit polyclonal antibody

别名： DBY

抗原： DDX3Y

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, and it has a homolog on the X chromosome. The gene mutation causes male infertility.

	ty, Sertoli cell-only syndrome or severe hypospermatogenesis, suggesting that this gene plays a key role in the spermatogenic process. Alternatively spliced variants, encoding the same protein, have been identified.
<b>Applications:</b>	IHC
<b>Name of antibody:</b>	DDX3Y
<b>Immunogen:</b>	Synthesized peptide derived from internal of human DDX3Y.
<b>Full name:</b>	DEAD (Asp-Glu-Ala-Asp) box helicase 3, Y-linked
<b>Synonyms :</b>	DBY
<b>SwissProt:</b>	O15523
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100

