

兔抗 CABIN1 多克隆抗体

中文名称：兔抗 CABIN1 多克隆抗体

英文名称： Anti-CABIN1 rabbit polyclonal antibody

别名： CAIN; PPP3IN; KB-318B8.7

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CABIN1

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Calcineurin plays an important role in the T-cell receptor-mediated signal transduction pathway. The protein encoded by this gene binds specifically to the activated form of calcineurin and inhibits calcineurin-mediated signal transduction. The encoded protein is found in the nucleus and contains a leucine zipper domain as well as several PEST motifs, sequences which confer targeted degradation to those proteins which contain them. Alternative splicing results in multiple transcript variants encoding two different isoforms.
Applications:	ELISA, IHC

Name of antibody:	CABIN1
Immunogen:	Synthetic peptide of human CABIN1
Full name:	calcineurin binding protein 1
Synonyms:	CAIN; PPP3IN; KB-318B8.7
SwissProt:	Q9Y6J0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	40-200

