

TPT1 抗原（重组蛋白）

中文名称：TPT1 抗原（重组蛋白）

英文名称：TPT1 Antigen (Recombinant Protein)

别名：HRF; p02; p23; TCTP

储存 冷冻（-20℃）

相关类别 抗原

概述

Full length fusion protein

技术规格

Full name:	tumor protein, translationally-controlled 1
Synonyms:	HRF; p02; p23; TCTP
Swissprot:	P13693
Gene Accession:	BC003352
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Involved in calcium binding and microtubule stabilization.