

兔抗 CABP5 多克隆抗体

中文名称：兔抗 CABP5 多克隆抗体

英文名称：Anti-CABP5 rabbit polyclonal antibody

别名：calcium binding protein 5; CABP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CABP5

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The product of this gene belongs to a subfamily of calcium binding proteins, which share similarity to calmodulin. Calcium binding proteins are an important component of calcium mediated cellular signal transduction. Expression of this gene is retina-specific. The mouse homolog of this protein has been shown to express in the inner nuclear layer of the retina, suggested its role in neuronal functioning. The specific function of this gene is unknown.
Applications:	ELISA, IHC
Name of antibody:	CABP5

Immunogen:	Full length fusion protein
Full name:	calcium binding protein 5
Synonyms:	CABP3
SwissProt:	Q9NP86
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	50-300



