

## 兔抗 DDX19A 多克隆抗体

中文名称：兔抗 DDX19A 多克隆抗体

英文名称：Anti-DDX19A rabbit polyclonal antibody

别名：DDX19L; DDX19-DDX19L

相关类别：一抗

储存：冷冻（-20℃）

抗原：DDX19A

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	ATP-dependent RNA helicase involved in mRNA export from the nucleus. Rather than unwinding RNA duplexes, DDX19 functions as a remodeler of ribonucleoprotein particles, whereby proteins bound to nuclear mRNA are dissociated and replaced by cytoplasmic mRNA binding proteins (By similarity).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DDX19A
<b>Immunogen:</b>	Fusion protein corresponding to a region derived from 279-478 amino acids of human DDX19A
<b>Full name:</b>	DEAD-box helicase 19A

<b>Synonyms:</b>	DDX19L; DDX19-DDX19L
<b>SwissProt:</b>	Q9NUU7
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC Positive control:</b>	Human liver cancer; Human esophagus cancer
<b>IHC Recommended dilution:</b>	20-100



