

兔抗 DDX27 多克隆抗体

中文名称：兔抗 DDX27 多克隆抗体

英文名称：Anti-DDX27 rabbit polyclonal antibody

别名：DEAD-box helicase 27 DRS1; RHLP; Drs1p; PP3241; HSPC259; dJ686N3.1

抗原：DDX27

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Fusion protein of human DDX27
Applications:	ELISA, WB, IHC
Name of antibody:	DDX27
Immunogen:	Fusion protein of human DDX27
Full name:	DEAD-box helicase 27
Synonyms:	DRS1; RHLP; Drs1p; PP3241; HSPC259; dJ686N3.1
SwissProt:	Q96GQ7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer;
IHC Recommend dilution:	20-100
WB Predicted band size:	90 kDa

WB Positive control:	HEPG2 , Hela and K562 cell lysates
WB Recommended dilution:	200-1000



