

TSSK3 抗原（重组蛋白）

中文名称：TSSK3 抗原（重组蛋白）

英文名称： TSSK3 Antigen (Recombinant Protein)

别名： testis specific serine kinase 3; TSK3; SPOGA3; STK22C; STK22D

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-268 amino acids of human TSSK3

技术规格

Full name:	testis specific serine kinase 3
Synonyms:	TSK3; SPOGA3; STK22C; STK22D
Swissprot:	Q96PN8
Gene Accession:	BC035354
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either germ cell differentiation or mature sperm function.