

兔抗 DDX39A 多克隆抗体

中文名称：兔抗 DDX39A 多克隆抗体

英文名称： Anti-DDX39A rabbit polyclonal antibody

别名： DExD-box helicase 39A; BAT1; DDXL; BAT1L; DDX39; URH49

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DDX39A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the DEAD box protein family. These proteins are characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD) and are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure, such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of the DEAD box protein family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene is thought to play a role in the prognosis of patients with gastrointestinal stromal tumors. A pseudogene of this gene is present on chromosome 13. Alternate splicing results

	in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known.
Applications:	ELISA, WB, IHC
Name of antibody:	DDX39A
Immunogen:	Fusion protein of human DDX39A
Full name:	DExD-box helicase 39A
Synonyms:	BAT1; DDXL; BAT1L; DDX39; URH49
SwissProt:	O00148
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human lung cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	49 kDa
WB Positive control:	231 and Jurkat cell lysates
WB Recommended dilution:	500-2000



