

## 兔抗 DDX39B 多克隆抗体

中文名称：兔抗 DDX39B 多克隆抗体

英文名称：Anti-DDX39B rabbit polyclonal antibody

别名：BAT1; UAP56; D6S81E

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DDX39B

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a member of the DEAD box family of RNA-dependent ATPases that mediate ATP hydrolysis during pre-mRNA splicing. The encoded protein is an essential splicing factor required for association of U2 small nuclear ribonucleoprotein with pre-mRNA, and it also plays an important role in mRNA export from the nucleus to the cytoplasm. This gene belongs to a cluster of genes localized in the vicinity of the genes encoding tumor necrosis factor alpha and tumor necrosis factor beta.

<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	DDX39B
<b>Immunogen:</b>	Synthetic peptide of humanDDX39B
<b>Full name:</b>	DEAD (Asp-Glu-Ala-Asp) box polypeptide 39B
<b>Synonyms :</b>	BAT1; UAP56; D6S81E
<b>SwissProt:</b>	Q13838
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	49 kDa
<b>WB Positive control:</b>	A431, SKOV3, 231, Jurkat, HepG2 and Hela cells
<b>WB Recommended dilution:</b>	500-2000



