

兔抗 CACNB3 多克隆抗体

中文名称：兔抗 CACNB3 多克隆抗体

英文名称： Anti-CACNB3 rabbit polyclonal antibody

别名： CAB3; CACNLB3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CACNB3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a regulatory beta subunit of the voltage-dependent calcium channel. Beta subunits are composed of five domains, which contribute to the regulation of surface expression and gating of calcium channels and may also play a role in the regulation of transcription factors and calcium transport.
Applications:	ELISA, IHC
Name of antibody:	CACNB3
Immunogen:	Synthetic peptide of human CACNB3
Full name:	calcium voltage-gated channel auxiliary subunit beta 3
Synonyms:	CAB3; CACNLB3

SwissProt:	P54284
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	40-200

