

TXN2 抗原（重组蛋白）

中文名称：TXN2 抗原（重组蛋白）

英文名称：TXN2 Antigen (Recombinant Protein)

别名：MTRX; TRX2; MT-TRX

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	thioredoxin 2
Synonyms:	MTRX; TRX2; MT-TRX
Swissprot:	Q99757
Gene Accession:	BC013726
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This nuclear gene encodes a mitochondrial member of the thioredoxin family, a group of small multifunctional redox-active proteins. The encoded protein may play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis.