

## 兔抗 CACNG1 多克隆抗体

中文名称：兔抗 CACNG1 多克隆抗体

英文名称：Anti-CACNG1 rabbit polyclonal antibody

别名：CACNLG

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CACNG1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

#### Background:

Voltage-dependent calcium channels are composed of five subunits. The protein encoded by this gene represents one of these subunits, gamma, and is one of two known gamma subunit proteins. This particular gamma subunit is part of skeletal muscle 1,4-dihydropyridine-sensitive calcium channels and is an integral membrane protein that plays a role in excitation-contraction coupling. This gene is part of a functionally diverse eight-member protein subfamily of the PMP-22/EMP/MP20 family and is located in a cluster with two family members that function as transmembrane AMPA receptor regulatory protein

	s (TARPs).
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CACNG1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human CACNG1.
<b>Full name:</b>	calcium channel, voltage-dependent, gamma subunit 1
<b>Synonyms :</b>	CACNLG
<b>SwissProt:</b>	Q06432
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	Hela cells lysate
<b>WB Recommended dilution:</b>	500-3000

