

TXNDC17 抗原（重组蛋白）

中文名称：TXNDC17 抗原（重组蛋白）

英文名称：TXNDC17 Antigen (Recombinant Protein)

别名：TRP14; TXNL5

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	thioredoxin domain containing 17
Synonyms:	TRP14; TXNL5
Swissprot:	Q9BRA2
Gene Accession:	BC006405
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Disulfide reductase. May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions. Modulates TNF-alpha signaling and NF-kappa-B activation. Has peroxidase activity and may contribute to the elimination of cellular hydrogen peroxide.