

## 兔抗 CACNG3 多克隆抗体

中文名称：兔抗 CACNG3 多克隆抗体

英文名称： Anti-CACNG3 rabbit polyclonal antibody

别名： calcium voltage-gated channel auxiliary subunit gamma 3;

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CACNG3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a type I transmembrane AMPA receptor regulatory protein (TARP). TARPs regulate both trafficking and channel gating of the AMPA receptors. This gene is part of a functionally diverse eight-member protein subfamily of the PMP-22/EMP/MP20 family. This gene is a susceptibility locus for childhood absence epilepsy.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	CACNG3
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	calcium voltage-gated channel auxiliary subunit gamma

	3
<b>SwissProt:</b>	O60359
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	Mouse kidney tissue lysate
<b>WB Recommended dilution:</b>	200-1000

