

TXNDC5 抗原（重组蛋白）

中文名称：TXNDC5 抗原（重组蛋白）

英文名称：TXNDC5 Antigen (Recombinant Protein)

别名：MUTED-TXNDC5

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	BLOC1S5-TXNDC5 readthrough (NMD candidate)
Synonyms:	MUTED-TXNDC5
Swissprot:	Q8NBS9
Gene Accession:	BC052310
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This locus represents naturally occurring read-through transcription between the neighboring MUTED (muted homolog) and TXNDC5 (thioredoxin domain containing 5) genes on chromosome 6. The read-through transcript is a candidate for nonsense-mediated mRNA decay (NMD) and is unlikely to produce a protein product.