

兔抗 CACNG5 多克隆抗体

中文名称：兔抗 CACNG5 多克隆抗体

英文名称： Anti-CACNG5 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： CACNG5

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a type II transmembrane AMPA receptor regulatory protein (TARP). TARPs regulate both trafficking and channel gating of the AMPA receptors. This gene is part of a functionally diverse eight-member protein subfamily of the PMP-22/EMP/MP20 family and is located in a cluster with two family members, a type I TARP and a calcium channel gamma subunit. This gene is a susceptibility locus for schizophrenia and bipolar disorder.
Applications:	ELISA, IHC
Name of antibody:	CACNG5
Immunogen:	Synthetic peptide of human CACNG5
Full name:	calcium voltage-gated channel auxiliary subunit gamma 5

SwissProt:	Q9UF02
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	40-200



