

TYROBP 抗原（重组蛋白）

中文名称：TYROBP 抗原（重组蛋白）

英文名称：TYROBP Antigen (Recombinant Protein)

别名：DAP12; KARAP; PLOSL

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	TYRO protein tyrosine kinase binding protein
Synonyms:	DAP12; KARAP; PLOSL
Swissprot:	O43914
Gene Accession:	BC011175
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a transmembrane signaling polypeptide which contains an immunoreceptor tyrosine-based activation motif (ITAM) in its cytoplasmic domain. The encoded protein may associate with the killer-cell inhibitory receptor (KIR) family of membrane glycoproteins and may act as an activating signal transduction element. This protein may bind zeta-chain (TCR) associated protein kinase 70kDa (ZAP-70) and spleen tyrosine kinase (SYK) and play a role in signal transduction, bo

ne modeling, brain myelination, and inflammation. Mutations within this gene have been associated with polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy (PLOS), also known as Nasu-Hakola disease. Its putative receptor, triggering receptor expressed on myeloid cells 2 (TREM2), also causes PLOS. Multiple alternative transcript variants encoding distinct isoforms have been identified for this gene.