

## 兔抗 DDX53 多克隆抗体

中文名称：兔抗 DDX53 多克隆抗体

英文名称： Anti-DDX53 rabbit polyclonal antibody

别名： CAGE; CT26

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： DDX53

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	DDX53
<b>Immunogen:</b>	Fusion protein of human DDX53
<b>Full name:</b>	DEAD (Asp-Glu-Ala-Asp) box polypeptide 53
<b>Synonyms :</b>	CAGE; CT26
<b>SwissProt:</b>	Q86TM3
<b>ELISA Recommended dilution:</b>	2000-10000

<b>IHC positive control:</b>	Human liver cancer and human gastric cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	71 kDa
<b>WB Positive control:</b>	HepG2 cells
<b>WB Recommended dilution:</b>	200-1000



