

UBE2S 抗原（重组蛋白）

中文名称：UBE2S 抗原（重组蛋白）

英文名 UBE2S Antigen (Recombinant Protein)

别名 EPF5; E2EPF; E2-EPF

储存 冷冻（-20℃）

相关类别 抗原

概述

Full length fusion protein

技术规格

Full name:	ubiquitin-conjugating enzyme E2S
Synonyms:	EPF5; E2EPF; E2-EPF
Swissprot:	Q16763
Gene Accession:	NP_055316
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the ubiquitin-conjugating enzyme family. The encoded protein is able to form a thiol ester linkage with ubiquitin in a ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins.