

UBE2T 抗原（重组蛋白）

中文名称：UBE2T 抗原（重组蛋白）

英文名称：UBE2T Antigen (Recombinant Protein)

别名：FANCT; PIG50; HSPC150

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	ubiquitin conjugating enzyme E2T
Synonyms:	FANCT; PIG50; HSPC150
Swissprot:	Q9NPD8
Gene Accession:	BC004152
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.