

兔抗 CAGE1 多克隆抗体

中文名称：兔抗 CAGE1 多克隆抗体

英文名称：Anti-CAGE1 rabbit polyclonal antibody

别名：cancer antigen 1; CT3; CT95; CTAG3; bA69L16.7

抗原：CAGE1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Cancer-associated gene 1 protein (CAGE1) is initially identified as a cancer/testis antigen. CAGE1 is mainly a nuclear protein that displays high expression in testis among normal tissues, and wide expression in various cancer tissues and cancer cell lines. CAGE1 may be associated with tumor progression, while the exact function of CAGE-1 is not well understood.
Applications:	ELISA, WB, IHC
Name of antibody:	CAGE1
Immunogen:	Fusion protein of human CAGE1
Full name:	cancer antigen 1

Synonyms:	CT3; CT95; CTAG3; bA69L16.7
SwissProt:	Q8TC20
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer;
IHC Recommend dilution:	25-100
WB Predicted band size:	90 kDa
WB Positive control:	Jurkat cell lysate
WB Recommended dilution:	200-1000



