

UBL3 抗原（重组蛋白）

中文名称：UBL3 抗原（重组蛋白）

英文名称：UBL3 Antigen (Recombinant Protein)

别名：HCG-1; PNSC1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	ubiquitin-like 3
Synonyms:	HCG-1; PNSC1
Swissprot:	O95164
Gene Accession:	BC059385
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Ubiquitin is a 77 amino acid protein that targets proteins for degradation by the 26S proteasome. Ubiquitin-like proteins are not directly involved in protein degradation, but appear to have many mechanistic similarities with the ubiquitin pathway. UBL3 (ubiquitin-like protein 3), also known as membrane-anchored ubiquitin-fold protein (MUB) or PNSC1, is a 117 amino acid membrane protein belonging to the ubiquitin-like family. Highly conserved between spec

ies, UBL3 contains two potential N-glycosylation sites, a potential p
rotein kinase C phosphorylation site, and a potential C-terminal pr
enylation site. The gene encoding UBL3 is localized to chromosom
e 13q12.3.