

兔抗 CALCOCO2 多克隆抗体

中文名称：兔抗 CALCOCO2 多克隆抗体

英文名称： Anti-CALCOCO2 rabbit polyclonal antibody

别名： calcium binding and coiled-coil domain 2; NDP52

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CALCOCO2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, WB
Name of antibody:	CALCOCO2
Immunogen:	Fusion protein of human CALCOCO2
Full name:	calcium binding and coiled-coil domain 2
Synonyms:	NDP52

SwissProt:	Q13137
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	52 kDa
WB Positive control:	Jurkat, Hela and A549 cell lysates
WB Recommended dilution:	200-1000

