

## 兔抗 CALHM1 多克隆抗体

中文名称：兔抗 CALHM1 多克隆抗体

英文名称： Anti-CALHM1 rabbit polyclonal antibody

别名： calcium homeostasis modulator 1; FAM26C

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CALHM1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a calcium channel that plays a role in processing of amyloid-beta precursor protein. A polymorphism at this locus has been reported to be associated with susceptibility to late-onset Alzheimer's disease in some populations, but the pathogenicity of this polymorphism is unclear.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CALHM1
<b>Immunogen:</b>	Fusion protein of human CALHM1
<b>Full name:</b>	calcium homeostasis modulator 1
<b>Synonyms:</b>	FAM26C

<b>SwissProt:</b>	Q8IU99
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	30-150

