

## UQCRH 抗原（重组蛋白）

中文名称：UQCRH 抗原（重组蛋白）

英文名称：UQCRH Antigen (Recombinant Protein)

别名：QCR6; UQCR8

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	ubiquinol-cytochrome c reductase hinge protein
<b>Synonyms:</b>	QCR6; UQCR8
<b>Swissprot:</b>	P07919
<b>Gene Accession:</b>	BC001426
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This protein may mediate formation of the c complex between cytochromes c and c1.