

USE1 抗原（重组蛋白）

中文名称：USE1 抗原（重组蛋白）

英文名称：USE1 Antigen (Recombinant Protein)

别名：D12; P31; SLT1; MDS032

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	unconventional SNARE in the ER 1
Synonyms:	D12; P31; SLT1; MDS032
Swissprot:	Q9NZ43
Gene Accession:	BC006005
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	SNARE that may be involved in targeting and fusion of Golgi-derived retrograde transport vesicles with the ER.