

兔抗 DECR1 多克隆抗体

中文名称：兔抗 DECR1 多克隆抗体

英文名称： Anti-DECR1 rabbit polyclonal antibody

别名： DECR; NADPH; SDR18C1

抗原： DECR1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

克隆类型： rabbit polyclonal

标记物： Unconjugate

技术规格

Background:	This gene encodes an accessory enzyme which participates in the beta-oxidation and metabolism of unsaturated fatty acyl-CoA esters. Auxiliary enzyme of beta-oxidation. It participates in the metabolism of unsaturated fatty acyl-CoA esters having double bonds in both even- and odd-numbered positions. Catalyzes the NADP-dependent reduction of 2,4-dienoyl-CoA to yield trans-3-enoyl-CoA.
Applications:	ELISA, WB, IHC
Name of antibody:	DECR1
Immunogen:	Synthetic peptide of human DECR1

Full name:	2,4-dienoyl CoA reductase 1, mitochondrial
Synonyms :	DECR; NADPH; SDR18C1
SwissProt:	Q16698
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	36 kDa
WB Positive control:	Human endometrial carcinoma and colon cancer tissue, mouse heart tissue and PC3 cells
WB Recommended dilution:	500-2000



