

UXT 抗原（重组蛋白）

中文名称：UXT 抗原（重组蛋白）

英文名称：UXT Antigen (Recombinant Protein)

别名：STAP1; ART-27

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	ubiquitously-expressed, prefoldin-like chaperone
Synonyms:	STAP1; ART-27
Swissprot:	Q9UBK9
Gene Accession:	BC000720
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.