

## 兔抗 CALU 多克隆抗体

中文名称：兔抗 CALU 多克隆抗体

英文名称： Anti-CALU rabbit polyclonal antibody

相关类别： 一抗

抗 原： CALU

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The product of this gene is a calcium-binding protein localized in the endoplasmic reticulum (ER) and it is involved in such ER functions as protein folding and sorting. This protein belongs to a family of multiple EF-hand proteins (CERC) that include reticulocalbin, ERC-55, and Calb45 and the product of this gene. Alternatively spliced transcript variants encoding different isoforms have been identified
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CALU
<b>Immunogen:</b>	Fusion protein of human CALU
<b>Full name:</b>	calumenin
<b>SwissProt:</b>	O43852

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	37 kDa
<b>WB Positive control:</b>	Human placenta tissue
<b>WB Recommended dilution:</b>	500-2000



