

兔抗 CAMK1G 多克隆抗体

中文名称：兔抗 CAMK1G 多克隆抗体

英文名称： Anti-CAMK1G rabbit polyclonal antibody

别名： VWS1; CLICK3; CLICKIII; dJ272L16.1

抗原： CAMK1G

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known.
Applications:	ELISA, IHC
Name of antibody:	CAMK1G
Immunogen:	Fusion protein of human CAMK1G
Full name:	calcium/calmodulin-dependent protein kinase IG
Synonyms :	VWS1; CLICK3; CLICKIII; dJ272L16.1
SwissProt:	Q96NX5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human brain

IHC Recommend dilution:

25-100

