

## 兔抗 DEFA4 多克隆抗体

中文名称：兔抗 DEFA4 多克隆抗体

英文名称： Anti-DEFA4 rabbit polyclonal antibody

别名： HP4; DEF4; HP-4; HNP-4

抗原： DEFA4

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Defensins are a family of microbicidal and cytotoxic peptides thought to be involved in host defense. They are abundant in the granules of neutrophils and also found in the epithelia of mucosal surfaces such as those of the intestine, respiratory tract, urinary tract, and vagina. Members of the defensin family are highly similar in protein sequence and distinguished by a conserved cysteine motif. Several alpha defensin genes are clustered on chromosome 8. This gene differs from other genes of this family by an extra 83-base segment that is apparently the result of a recent duplication.

	tion within the coding region. The protein encoded by this gene, defensin, alpha 4, is found in the neutrophils; it exhibits corticostatic activity and inhibits corticotropin stimulated corticosterone production.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DEFA4
<b>Immunogen:</b>	Synthetic peptide of human DEFA4
<b>Full name:</b>	defensin, alpha 4, corticostatin
<b>Synonyms :</b>	HP4; DEF4; HP-4; HNP-4
<b>SwissProt:</b>	P12838
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human brain
<b>IHC Recommend dilution:</b>	50-200

