

## 兔抗 DEFA6 多克隆抗体

中文名称：兔抗 DEFA6 多克隆抗体

英文名称： Anti-DEFA6 rabbit polyclonal antibody

别名： DEF6; HD-6

抗原： DEFA6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Defensins are a family of microbicidal and cytotoxic peptides thought to be involved in host defense. They are abundant in the granules of neutrophils and also found in the epithelia of mucosal surfaces such as those of the intestine, respiratory tract, urinary tract, and vagina. Members of the defensin family are highly similar in protein sequence and distinguished by a conserved cysteine motif. Several alpha defensin genes appear to be clustered on chromosome 8.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DEFA6

<b>Immunogen:</b>	Fusion protein of human DEFA6
<b>Full name:</b>	defensin, alpha 6, Paneth cell-specific
<b>Synonyms :</b>	DEF6; HD-6
<b>SwissProt:</b>	Q01524
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	50-200

