

VGLL2 抗原（重组蛋白）

中文名称：VGLL2 抗原（重组蛋白）

英文名称：VGLL2 Antigen (Recombinant Protein)

别名：VGL2; VITO1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	vestigial like family member 2
Synonyms:	VGL2; VITO1
Swissprot:	Q8N8G2
Gene Accession:	BC069316
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein with a transcriptional enhancer factor 1 (TEF-1) interaction domain. The encoded protein may act as a co-factor of TEF-1 regulated gene expression during skeletal muscle development. Alternatively spliced transcript variants encoding multiple isoforms have been observed.