

## 兔抗 DENND1B 多克隆抗体

中文名称：兔抗 DENND1B 多克隆抗体

英文名称： Anti-DENND1B rabbit polyclonal antibody

别名： DENN domain containing 1B; FAM31B; C1ORF18; C1orf218

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： DENND1B

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Clathrin (see MIM 118955)-mediated endocytosis is a major mechanism for internalization of proteins and lipids. Members of the connectin family, such as DENND1B, function as guanine nucleotide exchange factors (GEFs) for the early endosomal small GTPase RAB35 (MIM 604199) and bind to clathrin and clathrin adaptor protein-2 (AP2; see MIM 601024). Thus, connectins link RAB35 activation with the clathrin machinery.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	DENND1B

<b>Immunogen:</b>	Fusion protein of human DENND1B
<b>Full name:</b>	DENN domain containing 1B
<b>Synonyms:</b>	FAM31B; C1ORF18; C1orf218
<b>SwissProt:</b>	Q6P3S1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	86 kDa
<b>WB Positive control:</b>	Human fetal brain tissue and mouse brain tissue lysates
<b>WB Recommended dilution:</b>	500-2000

