

VIL1 抗原（重组蛋白）

中文名称：VIL1 抗原（重组蛋白）

英文名称：VIL1 Antigen (Recombinant Protein)

别名：VIL; D2S1471

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	villin 1
Synonyms:	VIL; D2S1471
Swissprot:	P09327
Gene Accession:	BC017303
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a family of calcium-regulated actin-binding proteins. This protein represents a dominant part of the brush border cytoskeleton which functions in the capping, severing, and bundling of actin filaments. Two mRNAs of 2.7 kb and 3.5 kb have been observed; they result from utilization of alternate poly-adenylation signals present in the terminal exon.