

VPS25 抗原（重组蛋白）

中文名称：VPS25 抗原（重组蛋白）

英文名称： VPS25 Antigen (Recombinant Protein)

别名： DERP9; EAP20; FAP20

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	vacuolar protein sorting 25 homolog
<b>Synonyms:</b>	DERP9; EAP20; FAP20
<b>Swissprot:</b>	Q9BRG1
<b>Gene Accession:</b>	BC006282
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II (ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A pseudogene of this gene is present on chromosome 1.