

兔抗 CAMSAP1 多克隆抗体

中文名称: 兔抗 CAMSAP1 多克隆抗体

英文名称: Anti-CAMSAP1 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: CAMSAP1

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Probable microtubule-binding protein that plays a role in the regulation of cell morphology and cytoskeletal organization. Through interaction with spectrin may regulate neurite outgrowth.
Applications:	ELISA, IHC
Name of antibody:	CAMSAP1
Immunogen:	Synthetic peptide of human CAMSAP1
Full name:	calmodulin regulated spectrin associated protein 1
SwissProt:	Q5T5Y3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

