

VTI1B 抗原（重组蛋白）

中文名称：VTI1B 抗原（重组蛋白）

英文名称：VTI1B Antigen (Recombinant Protein)

别名：VTI1; VTI2; VTI1L; v-SNARE; vti1-rp1; VTI1-LIKE

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	vesicle transport through interaction with t-SNAREs 1B
Synonyms:	VTI1; VTI2; VTI1L; v-SNARE; vti1-rp1; VTI1-LIKE
Swissprot:	Q9UEU0
Gene Accession:	BC003142
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	V-SNARE that mediates vesicle transport pathways through interactions with t-SNAREs on the target membrane. These interactions are proposed to mediate aspects of the specificity of vesicle trafficking and to promote fusion of the lipid bilayers. May be concerned with increased secretion of cytokines associated with cellular senescence.