

VWF 抗原（重组蛋白）

中文名称：VWF 抗原（重组蛋白）

英文名称：VWF Antigen (Recombinant Protein)

别名：VWD; F8VWF

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	von Willebrand factor
Synonyms:	VWD; F8VWF
Swissprot:	P04275
Gene Accession:	BC022258
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The glycoprotein encoded by this gene functions as both an antihemophilic factor carrier and a platelet-vessel wall mediator in the blood coagulation system. It is crucial to the hemostasis process. Mutations in this gene or deficiencies in this protein result in von Willebrand's disease. An unprocessed pseudogene has been found on chromosome 22.