

兔抗 CANT1 多克隆抗体

中文名称：兔抗 CANT1 多克隆抗体

英文名称： Anti-CANT1 rabbit polyclonal antibody

别名： calcium activated nucleotidase 1; DBQD; DBQD1; SCAN1; SHAPY; SCAN-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CANT1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This protein encoded by this gene belongs to the apyrase family. It functions as a calcium-dependent nucleotidase with a preference for UDP. Mutations in this gene are associated with Desbuquois dysplasia with hand anomalies. Alternatively spliced transcript variants have been noted for this gene.
Applications:	ELISA, WB
Name of antibody:	CANT1
Immunogen:	Fusion protein of human CANT1
Full name:	calcium activated nucleotidase 1
Synonyms:	DBQD; DBQD1; SCAN1; SHAPY; SCAN-1

SwissProt:	Q8WVQ1
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	45 kDa
WB Positive control:	Hela and A375 cell lysates
WB Recommended dilution:	500-2000

