

WBP1 抗原（重组蛋白）

中文名称：WBP1 抗原（重组蛋白）

英文名称：WBP1 Antigen (Recombinant Protein)

别名：WBP-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	WW domain binding protein 1
Synonyms:	WBP-1
Swissprot:	Q96G27
Gene Accession:	BC010012
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The globular WW domain, named for the conserved tryptophan residues in the protein motif present in various structural and regulatory proteins, is known to play a role in the mediation of protein-protein interactions. This gene encodes a ligand of the WW domain of the Yes kinase-associated protein. Readthrough transcription of the neighboring upstream gene, which encodes INO80 complex subunit B, into this gene generates a non-coding transcript.