

兔抗 DES 多克隆抗体

中文名称：兔抗 DES 多克隆抗体

英文名称：Anti-DES rabbit polyclonal antibody

别名：CSM1; CSM2; LGMD2R

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：DES

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a muscle-specific class III intermediate filament. Homopolymers of this protein form a stable intracytoplasmic filamentous network connecting myofibrils to each other and to the plasma membrane. Mutations in this gene are associated with desmin-related myopathy, a familial cardiac and skeletal myopathy (CSM), and with distal myopathies.
Applications:	ELISA, WB, IHC
Name of antibody:	DES
Immunogen:	Synthetic peptide of human DES
Full name:	desmin

Synonyms :	CSM1; CSM2; LGMD2R
SwissProt:	P17661
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	54 kDa
WB Positive control:	Jurkat and Hela cell
WB Recommended dilution:	200-1000



