

## 兔抗 CAP1 多克隆抗体

中文名称：兔抗 CAP1 多克隆抗体

英文名称： Anti-CAP1 rabbit polyclonal antibody

别名： CAP, CAP1-PEN

抗原： CAP1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is related to the <i>S. cerevisiae</i> CAP protein, which is involved in the cyclic AMP pathway. The human protein is able to interact with other molecules of the same protein, as well as with CAP2 and actin. Alternatively spliced transcript variants have been identified.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CAP1
<b>Immunogen:</b>	Fusion protein of human CAP1
<b>Full name:</b>	Adenylyl cyclase-associated protein 1

<b>Synonyms :</b>	CAP, CAP1-PEN
<b>SwissProt:</b>	Q01518
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human skin
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	Hela and HepG2 cells
<b>WB Recommended dilution:</b>	200-1000



