

WDR61 抗原（重组蛋白）

中文名称：WDR61 抗原（重组蛋白）

英文名称：WDR61 Antigen (Recombinant Protein)

别名：SKI8; REC14

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	WD repeat domain 61
Synonyms:	SKI8; REC14
Swissprot:	Q9GZS3
Gene Accession:	BC010080
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	WDR61 is a subunit of the human PAF and SKI complexes, which function in transcriptional regulation and are involved in events downstream of RNA synthesis, such as RNA surveillance (Zhu et al., 2005 [PubMed 16024656]).